

**Using Anesthesia Machines as ICU ventilators**  
**Specific to use at BIDMC**


**A. SETUP**

Insure manual ventilation device readily available – AMBU bag, PEEP valve, HMEF filter and Kelly Clamp

1. Connect/Check Central Gas Supplies
  - Check Line pressure – 45 psi or better o Full E-cylinders of oxygen and air as backup
  - Remove nitrous oxide hoses and cylinders
  - Vaporizers →Remove or drain (if not used for sedation)
  - Scavenging (White or purple outlet-depends on location)
    - o If not using the inhaled anesthetics – not necessary
    - o If performing volatile sedation – necessary
2. Configure machine with disposables Breathing Circuit o Filters
  - HMEF on airway, gas sampling on machine side
  - Second filter on the expiratory limb if possible.
  - Large (3 Liter) Reservoir Bag
  - Connect circuit to Gas analyzer for oxygen and carbon dioxide (use filter)
3. Perform Self-Test (q24hrs, but can be extended to q72hrs if needed)
  - Compliance measurement essential – do not change disposables after this
  - Confirm no errors
  - Check alarms, set limits, set to max volume (Defaults are not appropriate for ICU pts).
  - Set APL valve to 0 cmH2O

**Prevention of Contamination**

- Use both an HME (Y piece) and a viral filter (Expiratory limb)
- Gas sampling attached on the machine side of the HME filter
- Use of HME filter plus antiviral filter should be effective
- Should not need to decontaminate internal components



Preferred Filter Configuration  
 VFE > 99.99% for each filter. Gas sampling on machine side of filter. (Courtesy Draeger Medical)

**B. INITIATE THERAPY**

1. Fresh Gas Flow Options
  - Option 1: Low fresh gas flow to conserve oxygen (**not recommended for COVID pts → clogs filters**)
    - o Preserves humidity
    - o CO2 Absorbent must be available and maintained
    - o Inspired CO2 Alarm must be set to 5 mmHg
  - Option 2: Fresh gas flow => minute ventilation (1x – 1.5x MVV)
    - o No CO2 Absorbent needed (increase FGF if Inspired CO2 present)
2. Setting Oxygen Concentration on Electronic Flowmeters
  - i. Air/oxygen mix needed for delivered O2 concentration (see table)
3. Set Ventilator →Ventilation Mode & Settings
  1. Rate
  2. Volume
  3. I:E Ratio
  4. PEEP
4. Start Ventilator.

FIO2 Table.

|      | Oxygen to Air Ratio | Oxygen flow for 5 L/min total | Air flow for 5 L/min total |
|------|---------------------|-------------------------------|----------------------------|
| 21%  | 0 to 1              | 0.0                           | 5.0                        |
| 25%  | 0.06 to 1           | 0.3                           | 4.7                        |
| 30%  | 0.13 to 1           | 0.6                           | 4.4                        |
| 35%  | 0.21 to 1           | 0.9                           | 4.1                        |
| 40%  | 0.31 to 1           | 1.2                           | 3.8                        |
| 50%  | 0.59 to 1           | 1.9                           | 3.1                        |
| 60%  | 0.99 to 1           | 2.5                           | 2.5                        |
| 80%  | 3 to 1              | 3.8                           | 1.3                        |
| 100% | 1 to 0              | 5.0                           | 0.0                        |